

We leave the lifting to the lift. For residents who need assistance in standing, Roze stand-up model lifts are expressly designed to promote independence for the resident while at the same time relieving the lifting burden from the operator.



Strength Capacity

The unique height-adjustable kneepad and 450 pound weight capacity easily accommodates the majority of residents in your facility.



Ease of Use

With a large push handle, lighter steel and aluminum construction, and simple sling attachment, it is easy to accomplish everyday transfers.



Accommodating

For those who only need help rising, simply remove the footplate.



Maneuverability

A small footprint combined with lighter weight aluminum and steel construction provide easy maneuverability.



Invacare® Roze™ Stand-Up Lift

Height at sling hook-up - MAX	63.5"
Height at sling hook-up - MIN	37"
Base width OPEN	44"
Base width CLOSED	26"
Base height (Clearance)	4.8"
Base length	38.8"
Overall height	49.1"
Overall length	42"
Overall width	26"
Caster Size (Front/Rear)	4.0"
Weight capacity	450 lb.
Weight out-of-carton	117 lb.

Battery	24V DC (Rechargeable) 2 Shipped with Unit
Charger input	100-240V AC
Charger output/Charging time	29.5 DC Max 6 hours
Audio low battery alarm	Yes
Motor safety devices	Anti-Entrapment
Limited warranty electronics	3 years on Lift, 2 years on Electronics
Emergency stop button	Yes
Emergency up and down	Yes